

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

**Finland** 

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

## **Product Datasheet**

## HiFi DNA Polymerase Hot-Start #3000

Product Name HiFi DNA Polymerase Hot-Start

Catalog Number #3000

Description

HiFi DNA Polymerase Hot-Start is an aptamer-based hot-start, high fidelity hyperthermophilic recombinant DNA polymerase from the archaeon Pyrococcus furiosus. HiFi DNA Polymerase Hot-Start exhibits 5' to 3' polymerase activity and 3' to 5' exonuclease activity. The optimized buffer chemistry facilitates high sensitivity, yield, specificity, robust and rapid polymerase processivity. The enzyme is ideal for long, complex, difficult DNA templates and is resistant to PCR inhibitors. HiFi DNA Polymerase Hot-Start has a lower error rate than standard Taq DNA

polymerase (100x) and therefore, is suitable for PCR applications

where higher accuracy is needed.

Tested Applications

End-Point, Real-Time

Brand Medix Biochemica

Storage -20°C

Content 100 U: 2 U/µl, 1 x 0.05 ml HiFi DNA Polymerase Hot-Start; 1 x 1.7

ml 5x HiFi Buffer; 1 x 1.7 ml 10x HiFiOpt Enhancer

500 U: 2 U/μl, 1 x 0.25 ml HiFi DNA Polymerase Hot-Start; 3 x 1.7

ml 5x HiFi Buffer; 2 x 1.7 ml 10x HiFiOpt Enhancer