

Product Datasheet

Anti-h MPO 1701 SPTN-5 100266

| | |
|----------------------------|--|
| Product Name | Anti-h MPO 1701 SPTN-5 |
| Catalog Number | 100266 |
| Description | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Tested Applications | FIA |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/- 10 %) |
| Storage | +2-8°C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Epitope | N/D |
| Purity | ≥ 95 % |
| Affinity constant | KD= 1 x 10 ⁻⁹ M; KA= 1 x 10 ⁹ 1/M |
| Associated Products | Native MPO antigen Lee Biosolutions 426-10 |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| IEF Profile | 6.0-6.8 |
| Cross Reactivity | N/D |
| Specificity | Antibody recognizes human myeloperoxidase |
| Shelf Life | 36 months |