

Product Datasheet

Anti-Human CD154 Monoclonal Antibody, Azide Free Clone B-B29 854.290

Product Name Anti-Human CD154 Monoclonal Antibody, Azide Free Clone B-B29

Catalog Number 854.290

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications FC

Immunogen CD40L transfected 3T6 cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity Recognises the CD40L, TRAP, T-BAM antigen, a 28/30/33 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells