

Product Datasheet

Anti-Human CD45 Monoclonal Antibody, Azide Free Clone B-A11 854.300

Product Name Anti-Human CD45 Monoclonal Antibody, Azide Free Clone B-A11

Catalog Number 854.300

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications FC

Immunogen Thymus cells and Jurkat cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the Leucocyte Common Antigen (LCA) with different isoforms of 180-190-205-220 kDa

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells