

Product Datasheet

Anti-Human CD45 Monoclonal Antibody, Unconjugated Clone B-A11 854.303.020

Product Name	Anti-Human CD45 Monoclonal Antibody, Unconjugated Clone B-A11
Catalog Number	854.303.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	IHC, FC
Immunogen	Thymus cells and Jurkat cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the Leucocyte Common Antigen (LCA) with different isoforms of 180-190-205-220 kDa
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells