

Product Datasheet

Anti-h SAA 2201 SPTN-5 100279

| | |
|----------------------------|---|
| Product Name | Anti-h SAA 2201 SPTN-5 |
| Catalog Number | 100279 |
| Description | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Tested Applications | LF, IT, FIA |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/- 10 %) |
| Storage | +2-8°C |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| Epitope | Conformational |
| Purity | ≥ 95 % |
| Associated Products | Native human SAA, EastCoast Bio cat LA449, recombinant human SAA, Medix Biochemica cat 610070 and Lee Biosolutions cat 527-18 |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| IEF Profile | 7.4 - 8.8 |
| Cross Reactivity | Recognizes human SAA1, human SAA2, feline, canine and equine SAA. Does not recognize human SAA4 or bovine SAA. |
| Specificity | Antibody recognizes human serum amyloid A (SAA) |
| Shelf Life | 18 months |