

Product Datasheet

Anti-Human CD75 Monoclonal Antibody, Azide Free Clone B-L5 854.420

Product Name	Anti-Human CD75 Monoclonal Antibody, Azide Free Clone B-L5
Catalog Number	854.420
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC
Immunogen	Raji cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgM Kappa light chain
Clonality	Monoclonal
Source	Raji cell line
Specificity	Recognises lactosamines
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells