

Product Datasheet

Anti-Human CD86 Monoclonal Antibody, Azide Free Clone B-T7 854.450

Product Name	Anti-Human CD86 Monoclonal Antibody, Azide Free Clone B-T7
Catalog Number	854.450
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, Functional assay, FC
Immunogen	B70 transfected P815 cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the B7.2, B70 antigen, a 80 kDa protein
Biological Activity	Inhibits mixed lymphocyte culture
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells