

Product Datasheet

Anti-Human CD109 Monoclonal Antibody, Azide Free Clone B-E47 854.570

Product Name Anti-Human CD109 Monoclonal Antibody, Azide Free Clone B-E47

Catalog Number 854.570

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications FC

Immunogen Eahy 926 cell line

Alternative Names 854.570.000 854.570.005

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises a 170 kDa protein

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells