

Product Datasheet

Anti-Human CD184 Monoclonal Antibody, Azide Free Clone B-R24 854.600

Product Name Anti-Human CD184 Monoclonal Antibody, Azide Free Clone B-R24

Catalog Number 854.600

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications Functional assay, FC

Immunogen CXCR4 transfected CHO cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity Recognises the CXCR4 antigen, a 39.7 kDa protein

Biological Activity Inhibits HIV-1 infection

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells