

Product Datasheet

Anti-h IL-6 2703 SPTN-5 100328

Product Name	Anti-h IL-6 2703 SPTN-5
Catalog Number	100328
Description	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
Tested Applications	LF, FIA
Brand	Medix Biochemica
Form/Appearance	Liquid, may turn slightly opaque during storage
Concentration	5.0 mg/ml (+/- 10 %)
Storage	+2-8°C
Isotype	IgG1
Clonality	Monoclonal
Epitope	N/D
Purity	≥ 95 %
Affinity constant	KA= 1.3×10^{11} 1/M; KD= 7.7×10^{-11} M (=0.08 nM), BLI (Octet RED96e)
Associated Products	Recombinant IL-6 antigen 710015 , 710039 , 710038 and 341-36 .
Buffer	50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative
IEF Profile	6.3-7.6
Cross Reactivity	N/D
Specificity	Antibody recognizes human interleukin-6 (IL-6)
Shelf Life	24 months