

Product Datasheet

Anti-Human CD263 Monoclonal Antibody, Azide Free Clone B-D44 854.890

Product Name	Anti-Human CD263 Monoclonal Antibody, Azide Free Clone B-D44
Catalog Number	854.890
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	Functional assay, FC
Immunogen	Recombinant human TRAIL R3/Fc chimera
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises Dcr1 (CD263), a 75.90 kDa protein
Biological Activity	Inhibits TRAIL binding to TRAIL R3
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells