

Product Datasheet

Anti-Human CD218a Monoclonal Antibody, Azide Free Clone B-E43 854.900

Product Name	Anti-Human CD218a Monoclonal Antibody, Azide Free Clone B-E43
Catalog Number	854.900
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, Functional assay, FC
Immunogen	Recombinant human IL-18Ra/Fc
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the IL-18Ra, a 63 kDa protein. Cross-reacts with mouse CDw218a protein in cytometry.
Biological Activity	Inhibits human IFN γ secretion induced by IL-18 stimulated cells
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells