

Product Datasheet

Anti-Human CD218b Monoclonal Antibody, Azide Free Clone B-B46 854.910

Product Name Anti-Human CD218b Monoclonal Antibody, Azide Free Clone B-B46

Catalog Number 854.910

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, Functional assay

Immunogen Recombinant human IL-18Rb/Fc

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises the IL-18Rb, a 65.1 kDa protein

Biological Activity Inhibits human IFN γ secretion induced by IL-18 stimulated cells

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells