

Product Datasheet

Anti-Human Perforin Monoclonal Antibody, PerCP Conjugated Clone B-D48 854.958.010

Product Name	Anti-Human Perforin Monoclonal Antibody, PerCP Conjugated Clone B-D48
Catalog Number	854.958.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PerCP
Tested Applications	FC
Immunogen	Natural human Perforin
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 5% BSA and 0.1% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes human Perforin, a 70 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells