

Product Datasheet

Anti-Human CD8 Monoclonal Antibody, Azide Free Clone B-Z31 854.960

Product Name	Anti-Human CD8 Monoclonal Antibody, Azide Free Clone B-Z31
Catalog Number	854.960
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IHC, FC
Immunogen	Crude thymus membrane fraction
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the T-cell coreceptor, a 32-34 kDa protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells