

Product Datasheet

Anti-Human IL-1b Monoclonal Antibody, Azide Free Clone B-A15 855.010

Product Name Anti-Human IL-1b Monoclonal Antibody, Azide Free Clone B-A15

Catalog Number 855.010

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, ELISpot, Functional assay, FC

Immunogen Recombinant human IL-1b

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm.
Carrier and preservative free

Storage Store at +2-8°C for 12 months. DO NOT FREEZE.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Source Recombinant antibody from the original clone

Specificity Recognises both natural (pro and mature) and recombinant human IL-1b

Biological Activity Inhibits IL-1β induced proliferation on D10S murine cell line