

Product Datasheet

Anti-Human IL-4 Monoclonal Antibody, Azide Free Clone B-S4 855.030

Product Name	Anti-Human IL-4 Monoclonal Antibody, Azide Free Clone B-S4
Catalog Number	855.030
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	IHC
Immunogen	Recombinant human IL-4
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-4
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells