

Product Datasheet

Anti-Human IL-6 Monoclonal Antibody, Azide Free Clone B-E8 855.050

| | |
|----------------------------|--------------------------------------------------------------------------------------------|
| Product Name | Anti-Human IL-6 Monoclonal Antibody, Azide Free Clone B-E8 |
| Catalog Number | 855.050 |
| Description | Mouse monoclonal |
| Species/Host | Mouse |
| Conjugation | Azide free |
| Tested Applications | ELISA, ELISpot, Functional assay, FC |
| Immunogen | Recombinant human IL-6 |
| Alternative Names | 855.050.000, 855.050.005 |
| Brand | Diaclone |
| Form/Appearance | Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free |
| Storage | Stable at +2-8°C for 12 months. For longer storage freeze aliquots. |
| Isotype | Mouse IgG1 Kappa light chain |
| Clonality | Monoclonal |
| Specificity | Recognises both natural and recombinant human IL-6 |
| Biological Activity | Inhibits IL-6 induced proliferation on XG1 Cell line |
| Hybridoma | Myeloma X63/AG.8653 x Balb/c spleen cells |