

Product Datasheet

Anti-Human IL-13 Monoclonal Antibody, Azide Free Clone B-B13 855.140

Product Name Anti-Human IL-13 Monoclonal Antibody, Azide Free Clone B-B13

Catalog Number 855.140

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, Functional assay, FC

Immunogen Recombinant human IL-13

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human IL-13

Biological Activity Inhibits IL-13 induced proliferation on TF1 cell line

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells