

Product Datasheet

Anti-Human IL-15 Monoclonal Antibody, Azide Free Clone B-E29 855.220

Product Name	Anti-Human IL-15 Monoclonal Antibody, Azide Free Clone B-E29
Catalog Number	855.220
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, Functional assay, FC
Immunogen	Recombinant human IL-15
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-15
Biological Activity	Inhibits IL-15 induced proliferation on Kit 225 cell line
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells