

Product Datasheet

Anti-Human c-Myc Monoclonal Antibody, Unconjugated Clone 9E10 855.313.020

Product Name	Anti-Human c-Myc Monoclonal Antibody, Unconjugated Clone 9E10
Catalog Number	855.313.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	WB
Immunogen	Peptide corresponding to residues 408-439 in the human c-Myc protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Working dilution 1/50
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Peptide corresponding to residues 408-439 in the human c-Myc protein
Specificity	Recognises purified c-Myc protein
Hybridoma	Myeloma SP2/O x Balb/c splenocytes