

## Product Datasheet

### Anti-Human TNFa Monoclonal Antibody, Azide Free Clone B-F7 855.330

**Product Name** Anti-Human TNFa Monoclonal Antibody, Azide Free Clone B-F7

**Catalog Number** 855.330

**Description** Mouse monoclonal

**Species/Host** Mouse

**Conjugation** Azide Free

**Tested Applications** ELISA, ELISpot

**Immunogen** Recombinant human TNF-a

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. DO NOT FREEZE.

**Isotype** Mouse IgG1 Kappa light chain

**Clonality** Monoclonal

**Specificity** Recognises both natural and recombinant human TNF-a

**Hybridoma** Myeloma X63/AG.8653 x Balb/c spleen cells