

Product Datasheet

Anti-Human IL-17F Monoclonal Antibody, Azide Free Clone B-F60 855.360

Product Name Anti-Human IL-17F Monoclonal Antibody, Azide Free Clone B-F60

Catalog Number 855.360

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, ELISpot, Functional assay

Immunogen Recombinant human IL-17F

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human IL-17F

Biological Activity Neutralises the bioactivity of IL-17F on HT1080 cell line (measurement of IL-6 production)

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells