

Product Datasheet

Anti-Human MICA Monoclonal Antibody, Unconjugated Clone B-Z26 855.373.020

Product Name	Anti-Human MICA Monoclonal Antibody, Unconjugated Clone B-Z26
Catalog Number	855.373.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human MICA
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human MICA
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells