

Product Datasheet

Anti-Human CD62P Monoclonal Antibody, Azide Free Clone B-G43 855.390

Product Name Anti-Human CD62P Monoclonal Antibody, Azide Free Clone B-G43

Catalog Number 855.390

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA

Immunogen Recombinant human CD62P

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant CD62P

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells