

## Product Datasheet

### Anti-h Myoglobin 7004 SP-5 100354

<b>Product Name</b>	Anti-h Myoglobin 7004 SP-5
<b>Catalog Number</b>	100354
<b>Description</b>	Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components.
<b>Tested Applications</b>	LF, IT, FIA
<b>Brand</b>	Medix Biochemica
<b>Form/Appearance</b>	Liquid, may turn slightly opaque during storage
<b>Concentration</b>	5.0 mg/ml (+/- 10 %)
<b>Storage</b>	+2-8°C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Epitope</b>	N/D
<b>Purity</b>	≥ 95 %
<b>Affinity constant</b>	KA= 7 x 10 <sup>9</sup> 1/M
<b>Associated Products</b>	Recombinant Myoglobin antigen, Medix Biochemica 610030 and native Myoglobin antigen Lee Biosolutions 431-11.
<b>Buffer</b>	0.9 % NaCl, 0.095 % NaN <sub>3</sub> as a preservative
<b>IEF Profile</b>	5.4-6.1
<b>Cross Reactivity</b>	Hemoglobin < 0.005 %
<b>Specificity</b>	Antibody recognizes human myoglobin
<b>Shelf Life</b>	24 months