

Product Datasheet

Anti-Human CD62P Monoclonal Antibody, Azide Free Clone B-F46 855.400

Product Name Anti-Human CD62P Monoclonal Antibody, Azide Free Clone B-F46

Catalog Number 855.400

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, FC

Immunogen Recombinant human CD62P

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises CD62P, a 140 kDa single chain type I transmembrane glycoprotein present in secretory α-granules in platelets, Weibel-Palade bodies in endothelial cells & megakaryocytes

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells