

Product Datasheet

Anti-Human NKp46 Monoclonal Antibody, Azide Free Clone B-L46 855.490

Product Name Anti-Human NKp46 Monoclonal Antibody, Azide Free Clone B-L46

Catalog Number 855.490

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications Functional assay

Immunogen Recombinant human NKp46/Fc (CD335)

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity Recognises natural and recombinant NKp46

Biological Activity Measured by the ability of the coated antibody to induce degranulation of activated NK cells and to enhance IFN γ secretion by non-activated NK cells

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells