

## Product Datasheet

### Anti-Human NKp46 Monoclonal Antibody, Azide Free Clone B-L46 855.490

**Product Name** Anti-Human NKp46 Monoclonal Antibody, Azide Free Clone B-L46

**Catalog Number** 855.490

**Description** Mouse monoclonal

**Species/Host** Mouse

**Conjugation** Azide free

**Tested Applications** Functional assay

**Immunogen** Recombinant human NKp46/Fc (CD335)

**Brand** Diaclone

**Form/Appearance** Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

**Storage** Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

**Isotype** Mouse IgG2a Kappa light chain

**Clonality** Monoclonal

**Specificity** Recognises natural and recombinant NKp46

**Biological Activity** Measured by the ability of the coated antibody to induce degranulation of activated NK cells and to enhance IFNγ secretion by non-activated NK cells

**Hybridoma** Myeloma X63/AG.8653 x Balb/c node cells