

Product Datasheet

Anti-Human MICB Monoclonal Antibody, Azide Free Clone B-D54 855.520

Product Name	Anti-Human MICB Monoclonal Antibody, Azide Free Clone B-D54
Catalog Number	855.520
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC
Immunogen	Recombinant human MICB
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human MICB
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells