

Product Datasheet

Anti-Human Granulysin Monoclonal Antibody, Azide Free Clone B-L38 855.540

Product Name	Anti-Human Granulysin Monoclonal Antibody, Azide Free Clone B-L38
Catalog Number	855.540
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA, FC
Immunogen	Recombinant human Granulysin (NKG5 and LAG-2)
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human Granulysin
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells