

Product Datasheet

Anti-Human MCP-1 Monoclonal Antibody, Azide Free Clone B-K41 855.570

Product Name Anti-Human MCP-1 Monoclonal Antibody, Azide Free Clone B-K41

Catalog Number 855.570

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA

Immunogen Recombinant human MCP-1

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant MCP-1 (Monocyte Chemotactic Protein1), also called CCL2 or JE

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells