

## Product Datasheet

### Anti-Human Galectin-3 Monoclonal Antibody, Unconjugated Clone B-S43 855.613.020

<b>Product Name</b>	Anti-Human Galectin-3 Monoclonal Antibody, Unconjugated Clone B-S43
<b>Catalog Number</b>	855.613.020
<b>Description</b>	Mouse monoclonal
<b>Species/Host</b>	Mouse
<b>Conjugation</b>	Purified
<b>Tested Applications</b>	FC
<b>Immunogen</b>	Recombinant human Galectin-3
<b>Brand</b>	Diaclone
<b>Form/Appearance</b>	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
<b>Storage</b>	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
<b>Note</b>	Use 10 µl to label 10 <sup>6</sup> cells or 100 µl of whole blood
<b>Isotype</b>	Mouse IgG2a Kappa light chain
<b>Clonality</b>	Monoclonal
<b>Specificity</b>	Recognises both natural and recombinant protein
<b>Hybridoma</b>	Myeloma X63/AG8. x Balb/c node cells