

Product Datasheet

Anti-Human IL-24 Monoclonal Antibody, Azide Free Clone B-C60 855.630

Product Name	Anti-Human IL-24 Monoclonal Antibody, Azide Free Clone B-C60
Catalog Number	855.630
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human IL-24
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG2b Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant human IL-24 (interleukin-24) also called MDA-7 or Th2-specific cytokine FISP
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells