

Product Datasheet

Anti-Human IL-24 Monoclonal Antibody, Azide Free Clone B-C51 855.640

Product Name Anti-Human IL-24 Monoclonal Antibody, Azide Free Clone B-C51

Catalog Number 855.640

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA

Immunogen Recombinant human IL-24

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2a Kappa light chain

Clonality Monoclonal

Specificity Recognises both natural and recombinant human IL-24 (interleukin-24) also called MDA-7 or Th2-specific cytokine FISP

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells