

Product Datasheet

Anti-Human IL-35 Monoclonal Antibody, Azide Free Clone B-P35 855.660

Product Name	Anti-Human IL-35 Monoclonal Antibody, Azide Free Clone B-P35
Catalog Number	855.660
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	ELISA
Immunogen	EBI3 recombinant protein and IL-35/Fc recombinant protein
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises both natural and recombinant protein (no cross reaction with IL-12 and IL-27)
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells