

Product Datasheet

Anti-Human GM-CSF Monoclonal Antibody, Azide Free Clone B-Z13 855.690

Product Name	Anti-Human GM-CSF Monoclonal Antibody, Azide Free Clone B-Z13
Catalog Number	855.690
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human GM-CSF
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Store at +2-8°C for 12 months. DO NOT FREEZE.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises recombinant and natural GM-CSF (Granulocytemacrophage colony-stimulating factor)
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells