

Product Datasheet

Anti-Human GM-CSF Monoclonal Antibody, Azide Free Clone B-S37 855.700

Product Name Anti-Human GM-CSF Monoclonal Antibody, Azide Free Clone B-S37

Catalog Number 855.700

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA

Immunogen Recombinant human GM-CSF

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Store at +2-8°C for 12 months. DO NOT FREEZE.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises recombinant and natural GM-CSF (Granulocytemacrophage colony-stimulating factor)

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells