

Product Datasheet

Anti-Human NKp44 Monoclonal Antibody, Azide Free Clone B-N44 855.810

Product Name Anti-Human NKp44 Monoclonal Antibody, Azide Free Clone B-N44

Catalog Number 855.810

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications FC

Immunogen Recombinant human NKp44-FC

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile filtered through 0.22µm. Carrier and preservative free.

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognises natural and recombinant NKp44

Hybridoma Myeloma X63/AG.8653 x Balb/c node cells