

Product Datasheet

Anti-Human CD161 Monoclonal Antibody, Azide Free Clone B-D51 855.820

Product Name	Anti-Human CD161 Monoclonal Antibody, Azide Free Clone B-D51
Catalog Number	855.820
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide Free
Tested Applications	FC
Immunogen	Recombinant human KLRB1
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognises CD161
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells