

Product Datasheet

Anti-Human IL-32 Monoclonal Antibody, Azide Free Clone B-K42 855.900

Product Name	Anti-Human IL-32 Monoclonal Antibody, Azide Free Clone B-K42
Catalog Number	855.900
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Azide free
Tested Applications	ELISA
Immunogen	Recombinant human IL-32 gamma
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes both natural and recombinant human IL-32 (all isoforms)
Hybridoma	Myeloma X63/AG.8653 x Balb/c node cells