

Product Datasheet

Anti-Human CD196 Monoclonal Antibody, Unconjugated Clone B-R35 855.913.020

Product Name	Anti-Human CD196 Monoclonal Antibody, Unconjugated Clone B-R35
Catalog Number	855.913.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	huCCR6 NIH-3T3 transfected cells
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Specificity	CD196
Hybridoma	SP2/O x Balb/c B cells