

Product Datasheet

Anti-h Transferrin 3105 SPTN-5 100440

| | |
|----------------------------|--|
| Product Name | Anti-h Transferrin 3105 SPTN-5 |
| Catalog Number | 100440 |
| Description | Monoclonal mouse antibody, cultured in vitro under conditions free from animal-derived components. |
| Tested Applications | FIA |
| Brand | Medix Biochemica |
| Form/Appearance | Liquid, may turn slightly opaque during storage |
| Concentration | 5.0 mg/ml (+/- 10 %) |
| Storage | +2-8°C |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Epitope | N/D |
| Purity | ≥ 95 % |
| Affinity constant | KA= 5.3 x 10 ¹¹ 1/M; KD= 1.9 x 10 ⁻¹² M (=0.002 nM) |
| Associated Products | Apotransferrin (Sigma, T2252) |
| Buffer | 50 mM Na-citrate, pH 6.0, 0.9 % NaCl, 0.095 % NaN ₃ as a preservative |
| IEF Profile | 6.5-7.6 |
| Cross Reactivity | Does not recognize lactoferrin or hemoglobin. |
| Specificity | Antibody recognizes human transferrin |
| Shelf Life | 36 months |