

Product Datasheet

Anti-Human CD112 Monoclonal Antibody, Azide Free Clone B-C12 857.000

Product Name Anti-Human CD112 Monoclonal Antibody, Azide Free Clone B-C12

Catalog Number 857.000

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide free

Tested Applications ELISA, IP, FC

Immunogen KG-1 cell line

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. For longer storage freeze aliquots.

Isotype Mouse IgG2b Kappa light chain

Clonality Monoclonal

Specificity Recognises Nectin-2, the polio retrovirus receptor-2 (PRR2)

Hybridoma Myeloma X63/AG.8653 x Balb/c spleen cells