

Product Datasheet

Anti-Human Neurotensin R1 Monoclonal Antibody, Unconjugated Clone B-N6 857.043.020

Product Name	Anti-Human Neurotensin R1 Monoclonal Antibody, Unconjugated Clone B-N6
Catalog Number	857.043.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	NTR transfected CHO cell line
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE. DO NOT CONCENTRATE (risk of precipitation).
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgM Kappa light chain
Clonality	Monoclonal
Specificity	Recognises the membrane protein
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells