

Product Datasheet

Anti-Human DLL4 Monoclonal Antibody, Azide Free Clone B-D59 857.050

Product Name Anti-Human DLL4 Monoclonal Antibody, Azide Free Clone B-D59

Catalog Number 857.050

Description Mouse monoclonal

Species/Host Mouse

Conjugation Azide Free

Tested Applications FC

Immunogen Recombinant human DLL4

Brand Diaclone

Form/Appearance Phosphate-buffered saline. Sterile-filtered through 0.22 µm. Carrier and preservative free

Storage Stable at +2-8°C for 12 months. DO NOT FREEZE.

Isotype Mouse IgG1 Kappa light chain

Clonality Monoclonal

Specificity Recognizes natural and recombinant human DLL4

Hybridoma Myeloma X63/AG.8653 x Balb/c lymph node cells