

Product Datasheet

Anti-Human DLL4 Monoclonal Antibody, Unconjugated Clone B-D59 857.053.020

Product Name	Anti-Human DLL4 Monoclonal Antibody, Unconjugated Clone B-D59
Catalog Number	857.053.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Recombinant human DLL4
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Specificity	Recognizes natural and recombinant human DLL4
Hybridoma	Myeloma X63/AG.8653 x Balb/c lymph node cells