

Product Datasheet

IgG1 Isotype Control, Unconjugated Clone B-Z1 857.073.020

Product Name	IgG1 Isotype Control, Unconjugated Clone B-Z1
Catalog Number	857.073.020
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	Purified
Tested Applications	FC
Immunogen	Rat monocytes and macrophages
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1% BSA and 0.09% sodium azide
Storage	Stable at +2-8°C for 12 months. For longer storage freeze aliquots.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG1 Kappa light chain
Clonality	Monoclonal
Source	Rat monocytes and macrophages
Specificity	Negative control for biological test
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells