

Product Datasheet

IgG2a Isotype Control, PerCP Conjugated Clone B-Z2 857.088.010

Product Name	IgG2a Isotype Control, PerCP Conjugated Clone B-Z2
Catalog Number	857.088.010
Description	Mouse monoclonal
Species/Host	Mouse
Conjugation	PercP
Tested Applications	FC
Immunogen	Rat monocytes and macrophages
Brand	Diaclone
Form/Appearance	Phosphate-buffered saline with 1 % BSA and 0.1 % sodium azide
Storage	Stable at +2-8°C for 12 months. DO NOT FREEZE.
Note	Use 10 µl to label 10 ⁶ cells or 100 µl of whole blood
Isotype	Mouse IgG2a Kappa light chain
Clonality	Monoclonal
Source	Rat monocytes and macrophages
Specificity	Negative control for biological test
Hybridoma	Myeloma X63/AG.8653 x Balb/c spleen cells